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1 - GENERAL INFORMATION

1.1 - CE mark

These compressors bear CE conformity mark according to directive 2006/42/EC, ensuring that they are designed and manufactured in compliance with the safety specifications of the European Community directives.

1.2 - General instruction

This manual is a part of the compressor and must always stored close to it even in case of its selling.

Compressor owner and/or operator must know the content of this manual in order to perform its correct use and its safe maintenance.

When operator is not able to totally understand this manual instructions, due to language reason, dealer must supply a complete and correct translation on operator first language.

The respect of the instructions contained on this manual ensures the unit safe use and its correct and reliable working.

The compressor installation must be performed only by skilled operators. Same for the maintenance tasks: in case of parts replacement, it can be used only original spare parts or parts having features equivalent to the original ones.

No modifications to the compressor configuration is approved; no setting changes is approved, except for the ones reported on the present manual.

The manufacturer shall not be held liable for any damage to person or objects due to an improper or not-permitted use of the compressor, or to lack or uncorrect maintenance performing.

1.3 - Safety instruction.

This compressor has been designed and manufactured in a way to avoid any risks when operated according to the present instructions; however the following warnings are reported to avoid any residual dangers:

- do not use the compressor for any other purposes than those it has been designed to;
- keep it in a covered place and protected from rain and humidity.
- install it out of reach of children and of un-skilled people;
- do not compress any other fluid than air;
- never use the compressor in an environment with an explosive atmosphere;

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- never repair the compressor while it is powered or when its tank is under pressure;
- always use genuine spare parts and components suited to the performances of this unit;
- never power the compressor by an other voltage than what indicated on the product plate;
- never power the compressor when the supply cable is damaged;
- never switch on the compressor when it looks damaged; stop operating the compressor when it's making strange noise.
- never carry the compressor without its original packaging.

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1.4 - Danger instructions and markings.

The following markings are used on the compressor and on this manual in order to remark sensitive instructions or situations:

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
Warning/Ban to prevent any damages to people or things

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
Warning for danger due to hot surface

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
Warning for automatic restarting of the unit

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5 - Product data plate.

It follows the structure of the data plate stucked on each product.

1	 Werther International S.p.A. via Brunelleschi 12 42124 Cadè Reggio Emilia - ITALY			5
			G/2010	6
2	1/24/3	V Hz ~ 220 60 1	 l/min cfm 100 3.57	7
3	Serial N. 000000056	W A 800 7.8	 l gal 24 6.34	8
4	 dB A 70	IP20	 bar psi 7 101	9
			MADE IN ITALY	10
	11	12	13	

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2 - PRODUCT INFORMATION

2.1. - Description

Compressor has been designed for the production of compressed air, to be used by other device.

It is mainly composed by an oli-less motor-compressor, a tank for the compressed air storage, a safety valve able to discharge the over pressure and a pressure-switch, for the automatic operating of the device.

Some models are equipped with a drying system described at point 2.4.

When compressor is employed to feed up medical devices, installation must be performed complying with the specifications of European Community Directive 94/32/EEC.

The compressed air produced by this unit is not suitable to feed directly any device having the purpose to use the air on the patient body or through any part in

physical contact with the patient body or for the human breathing.


Avoiding the safety indications listed at point 1-4 of the present manual results in an improper use and it can cause a highly risk of danger.


Avoiding this manual reading can result in a danger to people and things.

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The code that distinguishes every compressor model consists of 5 dials: **A/BB/CD**; below the meaning of the dials

A Motor type 1, 2, 3, 4

BB Tank capacity in litres: 24/50

C Tank form: 3 = horizontal
5 = pancake

D Absente = without dryer
7 = with dryer
9 = with sound-cabinet

Examples:

2/24/37 _____ with dryer
 _____ tank form: horizontal
 _____ tank capacity: 24 Lt.
 _____ motor type: 2

3/50/5 _____ without dryer
 _____ tank form: pancake
 _____ tank capacity: 50 Lt.
 _____ motor type: 3

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
2.4 - ADSORPTION DRYER

The equipment is composed by:

- Radiator with fan
- Drying column with desiccant material
- Outlet filters
- Electrovalve for column draining
- Tank for drainage
- Inlet filter

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DRYING PROCESS

The dryer is composed by three parts: inlet connection with a 5 micron filter, aluminium column filled with a bed of desiccant material, outlet with 5 micron filter.

In the bottom side of the column, an electrovalve operates the water draining collected during the regeneration process.

The working principle is based on the adsorption effect performed by the bed of desiccant material while it is crossed by the air flow. During its flowing from the bottom side of the column to the upper, the compressed air progressively loses its moisture, captured by the desiccant material, and the air humidity rate decreases.

When the motor stops running, the compressed air stops flowing and the column remains under pressure. At the same time, the electrovalve opens and the compressed air contained in the column flows downwards. So the desiccant bed is crossed by the air flow at a lower humidity rate and frees the captured moisture. Then the humid air flow is released outside through the electrovalve and the desiccant material recovers its adsorption feature.

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2.5 - PACKING

The compressor is shipped in a single pack made up of a pallet and a cardboard box.
Gross weight and packing dimensions of each model is reported on table at chapter 2.3.

2.6 - LIFTING AND HANDLING

The compressors have to be handled and positioned using fork-lift trucks or transpallets.

2.7 - STORAGE


The compressors packed have to be kept in a dry, covered and sheltered place at a temperature between -10°C and $+50^{\circ}\text{C}$ and a relative humidity between $+10\%$ and $+90\%$, with no condensation.

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Packages can be stacked but not more than 2 pieces.

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3 - INSTALLATION

3.1. - GENERAL INSTRUCTIONS

The compressor installation must be performed only by trained or skilled operators.

Take the compressor out of its package and check that no damages occurred during transport; then assemble the provided silent-feet (fig.1).

Install the compressor in a suitably sized room, well ventilated and not wet /relative humidity less than 75%), where the temperature is between +5°C and +40°C.

Install the compressor on a flat surface, made of a material able to support this device weight.

When this compressor is not placed on a floor, check that the frame supporting the compressor is properly sized to carry on the compressor weight and is steady to the dynamic stresses caused by its running.

In case of not enough ventilation, install a suitably sized exhaustor or fan.

Avoid to place the compressor close to walls or furnitures; leave a free space (almost 50 cm.) around the unit in order to ensure an adequate ventilation.

3.2 - ELECTRICAL WIRING

Connect the unit to an electrical line equipped with a ground cord and having a voltage matching the data reported on the product data plate.

A tolerance of +/-10% is allowed on the line voltage.

The electrical line must be protected by a device/circuit breaker properly sized according to the maximum power/current required by the compressor.

This protection device must be able to avoid breaking the circuit while the motor is starting and its current rating can be 5-6 times greater than the standard one.

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The following table shows the current rate for each motor type:

	115 V - 60Hz - 1~	220V - 60Hz - 1~	230V - 50Hz - 1~	400 V - 50Hz - 3~
Motore tipo 1 Motor type 1	800W - 7,8A	800W - 4,8A	1000W - 5,8A	-
Motore tipo 2 Motor type 2	1450W - 17A	1450W - 8,5A	1300W - 7,5A	-
Motore tipo 3 Motor type 3	1970W - 18A	1970W - 9A	2100W - 10A	2100W - 6A
Motore tipo 4 Motor type 4	-	2150W - 11A	2480W - 12A	2480W - 5,2A

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The power supply plug must comply with the standard in force in the Country where the compressor is put in use and it must be proper to the compressor maximum current rating.

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Check that the supplied electrical cable is not damaged; in case of cable extension, use flexible cable with earth cord type H05VV or SJT (or having better features), and sections according to the following table:

	Fino a 3 m. To 3 mt.	
	220 V 50 Hz	115 V 60 Hz
Motore tipo 1 Motor type 1	1 mm ²	1,5 mm ²
Motore tipo 2 Motor type 2	1,5 mm ²	2,5 mm ²
Motore tipo 3 Motor type 3	1,5 mm ²	2,5 mm ²
Motore tipo 4 Motor type 4	1,5 mm ²	2,5 mm ²


	da 3 a 20 m. From 3 to 20 mt.	
	220 V 50 Hz	115 V 60 Hz
Motore tipo 1 Motor type 1	1,5 mm ²	2,5 mm ²
Motore tipo 2 Motor type 2	2,5 mm ²	4 mm ²
Motore tipo 3 Motor type 3	2,5 mm ²	6 mm ²
Motore tipo 4 Motor type 4	2,5 mm ²	6 mm ²

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In case of triphase compressor, at the first start check the real direction of rotation of the cooling fan (see arrow on the casing).

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
Before the compressor is put in use, check that the ground circuit is properly connected.

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




3.4 - PNEUMATICAL CONNECTION

Perform the pneumatical connection to the compressor by a properly sized hose, according to the compressor air flow rate and its maximum pressure.

Perform the connection to the compressor in a safe and durable way by pneumatic connector, equipped with a locking system by a joint or by a thread and proper to the maximum pressure of the compressor.

Check the pneumatic line connected to the compressor is free of any air leakage; this event can result in a unexpected compressor starting.

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4 - OPERATING INSTRUCTIONS

Check that the mains voltage corresponds to the one showed on the product data plate.

Turn the pressure switch button to position "0" (fig.2).

Insert the plug into the outlet and start the compressor turning the pressure switch button to position "1".

The compressor working process is automatic.

The pressure switch stops the compressor when the pressure in the tank reaches the maximum value allowed and starts it again when the pressure drop to the minimum value.

Adjust the outlet air pressure operating on the reducer-filter button (fig.3):

- a) pull up the button;
- b) turning the button clockwise the pressure rise;
- c) turning the button counterclockwise the pressure decrease;
- d) once reached the suitable pressure, push the button down.

In case of compressor with dryer, check that the compressor water collecting hose is connected to the tank.

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Be aware that compressor motor head and its outlet hose can get high temperature. Do not touch these parts to avoid the risk of thermal shock.

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5 - MAINTENANCE

EVERY WEEK

- Drain the water collected in the air outlet filter (Fig. 4) proceeding as follows: (do these operations with the tank under pressure)
 - 1) turn the button clockwise;
 - 2) pushing the button drain the dew contained in the bowl;
 - 3) at the end of this operation turn the button counterclockwise.
- Drain trough the drain plug (Fig. 6) the dew contained in the tank and empty the tank.

EVERY MONTH

- Disassemble the air intake filter (Fig. 5) and clean it blowing in compressed air.
- Check the compressor efficiency; perform the following checks:
 - mechanical: check the screws tightening, the motor shock absorbers integrity, the lack of rust on the frame and on the tank; the efficiency of the motor ventilation system;
 - pneumatical: check the connectors tightening, the lack of damages on the pressurized hose, the integrity of the safety valve and its correct opening (remove any block to its proper opening), the lack of air leakage;
 - electrical: check the integrity of the electrical cord and its connections, the correct operating of the pressure switch, the efficiency of the ground protective circuit.Check the unit stability and verify the integrity of its system of fixing or positioning.

EVERY SIX MONTHS

- Disassemble the coalescing cartridge contained in the filter (Fig. 4) and clean it by blowing compressed air. replace the filter if needed. This operation must be done with tank completely out of pressure.
- Check the dryer performance and, if needed, replace the desiccant material cartridge.

EVERY YEAR

- Check the safety valve. Check the valve is not damaged and perform the control of its opening: pull the valve ring (Fig.7) and check the air releasing. Check the valve gets closed after the ring is set free.

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Be aware that compressor motor head and its outlet hose can get high temperature. Do not touch to avoid the risk of thermal shock.

The control and maintenance operations must be performed only by trained or skilled operators.

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6 - Troubles - Reasons - Remedies

THE COMPRESSOR STOPS WHILE WORKING

- Models with type 1 motor:

- This motor has an automatic resetting thermal protection, that stops the compressor when the temperature gets too high. The compressor will start again automatically after 15/20 minutes.

- Models with type 2,3 and 4 motor:

- These motors have a manual resetting thermal protection placed on the terminal box cover (Fig. 8). If while working the thermal protection starts operating due to overheating, wait for few minutes, then manually reset the thermal switch.

AIR LEAKS

- Can be due to any connection bad seal, check all connections wetting with suds.

THE COMPRESSOR DOES NOT REACH THE MAXIMUM PRESSURE.

- Check any air leak.
- Check the pressure switch efficiency and if necessary adjust it (Fig. 9). Remove the cover. To adjust the maximum pressure operate on screws "a" and "b". Turning clockwise the pressure rise, turning counterclockwise the pressure decrease.

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LEAK FROM THE VALVE PLACED UNDER THE PRESSURE SWITCH

- Damaged valve, replace it. The defect can be caused by the unperfect valve seal (Fig. 10). Empty completely the tank of the air contained, unscrew the valve cap, clean the rubber disk and his place with a dry cloth and assembly the whole with care.

THE COMPRESSOR RUNS BUT DOES NOT LOAD

- The defect can be due to the valve or the gasket break (Fig. 11). Replace immediately the damaged part.
- The defect can be due to the piston seal ring wear (Fig. 12). Replace immediately the damaged part.

THE COMPRESSOR CONTINUE TO DRAIN AIR THROUGH THE DRYER ELECTROVALVE

- Disassemble the electrovalve and blow in compressed air (Fig. 13); if necessary replace the electrovalve.

CONDANSATE INSIDE THE AIR RECEIVER (ONLY FOR MODEL EQUIPPED WITH DRYER)

- Replace the dryer cartridge

All the operation mentioned in chapter "TROUBLES - REASONS - REMEDIES" **must be done with the tank completely empty and with the plug disconnected.** It is necessary that all the operations are carried out by authorized people.

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7 - SPECIAL INFORMATION

7.1 - COMPRESSOR DEMOLITION

During the compressor demolition all possible safety regulations must be observed in order to avoid any damage to people or things.

All the metal parts can be recycled; rubber and plastic parts have to be disposed of in the compliance with the laws in force in the country where the compressor is installed.

The packing material has to be disposed of in compliance with the regulation in force in the country where the compressor is being installed, or recycled or reused.

7.2 - SPARES

Defective parts must be replaced only by people authorized by the retailer shown on the cover; all possible safety regulations must be observed in order to avoid any damage to people or things.

HOW TO ORDER SPARE PARTS

To order the spare parts it is necessary to: indicate the compressor serial number and the year of manufacturing as shown in the cover.

indicate the part code requested referring to exploded drawings for your type of compressor.

Among the codes you will find some of them ending by the letters XX (ex.: C0545XX).

These letters refer to the colour numbers listed here below. When you have to order a piece in your own colour, you simply have to replace the two XX with the colour number. indicate the quantity requested.

The request must be placed with the **authorized distributor** shown in the cover.

Colours table	
Code	Colour
01	Black
07	Gray RAL 7016
09	White RAL 9002

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














Failure to observe the safety regulations may cause serious damage to people or things. The manufacturer shall not be held liable for damage to people or things caused by improper or not permitted use of the compressor.


7.3 - WARRANTY





- The compressor is guarantee for 12 months from date of sale.
- It regards only the free replacement of parts recognised as defective by the manufacturer apart from the electric parts and worn parts.
- The guarantee automatically cease in case of tampering and bad usage.
- The warranty does not include transport and labour costs.

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